


Aquatec® | **Bathing Transfer aids**



Aquatec **Trans**

- Compatible with the *Aquatec* range of bath lifters
- Large comfortable seating surface
- Safe: turntable is only unlocked by the weight of the user
- Seat surface turns and slides easily
- Soft, skin-friendly cover: easy-care, wipes clean

3 year warranty



Aquatec **Disk**

- Large comfortable seating surface
- Turns very easily
- Soft, skin-friendly cover: easy-care, wipes clean
- The *Aquatec Disk* is designed to be used in conjunction with the *Aquatec Dot* shower stool

3 year warranty

AQUATEC®



Yes, you can.®



Aquatec **Step**

- Stable frame
- Feet with non-slip rubber caps
- Non-slip, removable cover, machine-washable

3 year warranty

150kg/24 stone weight limit

Technical data

Order numbers	Aquatec Trans	Aquatec Disk	Aquatec Step
Turning and transfer aid, blue	4.03.002		
Turning and transfer aid, grey	4.03.008		
Rotating disk, blue		4.02.002	
Rotating disk, grey		4.02.008	
Step, blue			3.00.002
Dimensions			
Diameter	365 mm	365 mm	–
Length	530 mm	–	400 mm
Width	–	–	300 mm
Height	25 mm	25 mm	165 mm
More info			
Weight	2.2 kg	1.6 kg	1.5 kg
Load capacity	130 kg (21 stone)	130 kg (21 stone)	150 kg (24 stone)
Colours	Blue, grey	Blue, grey	Blue, white
Warranty	3 Years	3 Years	3 Years