

KONCENTRATOR KISIKA ZA KUĆNU UPOTREBU

8F-10 JEDNOSMJERNI PROTOK



Medicinska kvaliteta



Visoka koncentracija kisika



Protok kisika do 10 L



Brojač sati rada



Detekcija kvara



Kotačići



HEPA filter

- Maksimalni preporučeni protok: 0-10 L/min
- Koncentracija kisika: 87 % - 96 %
- Maksimalni izlazni tlak: 60kPA
- Razina buke: 49 dB(A)
- Ulazna snaga: 700 VA
- Minimalno vrijeme rada: 15 minuta



CE